



FY 2026 ICD-10 Codes Installation Guide

Installing FY2026 ICD-10 Code Updates to the WinOMS Software

This job aid contains instructions for updating the WinOMS software with the International Classification of Diseases, Tenth Revision, Clinical Modification (ICD-10 CM) diagnosis code set additions for Fiscal Year 2026, using the ICD-10 Code Update Utility.

Note: The updated code set will be included in WinOMS, v9.8.x and v10.1.x when they become available. Once you install one of these versions, the utility is not needed.

To use this utility, you must:

- Be using WinOMS, version 8.11 or higher.
- Download and run this utility on a server computer.
- Have Administrator security permissions to install this update.

To download and run the ICD-10 Update Utility:

- 1 Close WinOMS on all machines.
- 2 On the server computer, go to <https://gosensei.com/cdt-icd-codes>.
- 3 Scroll down to **WinOMS ICD-10 Code Updates**, select your version: **v9** or **v10**, and click **Download Update**. A short form is displayed.
- 4 Complete the form and click **Submit**. The appropriate file is downloaded to your default download location, or you are prompted to save it, depending on your browser.
(Version 10.x: **V10_ICD10Update_26_01.7z**, Version 9.x: **ICD10Update_26_01.7z**)
- 5 Navigate to the folder and verify the correct file has been downloaded.
- 6 Right-click the file and select **Run as administrator**. The **ICD-10 Code Update Utility Setup Wizard** is displayed.
- 7 Click **Next**. The **Select Destination Location** window is displayed. The server location should be populated, but if not, click **Browse** and select **c:\WinOMSCS**.
- 8 Click **Next**. The **Ready to Install** window is displayed.
- 9 Click **Install**. The utility is installed, and you are prompted for login credentials.
- 10 Enter your WinOMS user ID and password, and then click **OK**. The **ICD-10 Code Update Utility FY2026** utility begins processing the code changes.
When the status reads **Process Complete**, ICD-10 codes have been updated in the database.

To confirm ICD-10 utility has executed the update:

Navigate to **c:\winomscs** and open the file: **ICD10Utility.log**.

To add new codes to your tables:

- 1 Select **Tables > Diagnosis > Codes > New**.
- 2 Under **Actions**, enter the code you want (i.e., **J34.82**) and click **Lookup**.
- 3 Select a **Category** and click **Save**. The code is added to the list.

If you have issues during the installation, contact a Carestream Dental Support representative at: <https://gosensei.com/winoms-support/>.